

Olerup SSP[®] HLA-B*73 Kit



The Olerup SSP[®] product line expands to include a 3 well PCR-SSP high resolution kit for the HLA-B*73 group.

- Ready-to-test, pre-loaded PCR wells complete with primers and all reagents
- Single-test strips for easy use

The newly assigned allele HLA-B*73:02 was described in the IMGT/HLA Database release 3.1.0, 16 July 2010¹. The HLA-B*73 group now consist of two alleles, the B*73:01 and the B*73:02.

Olerup SSP[®] offers a complete product line for high resolution HLA typing and the product line is continuously developed to meet the needs of the market.

Olerup SSP[®] is now introducing a new high resolution kit for the HLA-B*73 group. The kit consists of three wells; One HLA-B*73 group specific well and one well each for the B*73:01 and the B*73:02 alleles.

Ordering information:

Product	Contents	Product No. Including <i>Taq</i> Polymerase	Product No. Without <i>Taq</i> Polymerase
HLA-B*73	6 single-strip tests, Master Mix	101.555-06	101.555-06u

Ref 1. Robinson J, Waller MJ, Parham P, de Groot N, Bontrop R, Kennedy LJ, Stoehr P, Marsh SGE; IMGT/HLA and IMGT/MHC: sequence databases for the study of the major histocompatibility complex; Nucleic Acids Research (2003), 31:311-314

Manufactured by:
OLERUP SSP AB
Address: Hasselstigen 1
SE-133 33 Saltsjöbaden, Sweden
E-mail: info-ssp@olerup.com
Tel: +46 8 717 88 27. Fax: +46 8 717 88 18



www.olerup-ssp.com



Distributed by:

MILAN ANALYTICA AG
YOUR BIOTECHNOLOGY PARTNER

Baslerstrasse 15
4310 Rheinfelden
Switzerland
www.milananalytica.ch
info@milananalytica.ch

Tel.: +41 61 845 99 88
Fax: +41 61 845 99 87

Olerup SSP® AB is ISO 9001:2008 and ISO 13485:2003 certified.
Olerup SSP® AB is IVD certified in accordance with the requirements of the IVD Directive 98/79/EC.
Olerup SSP® HLA Typing Kits are FDA cleared and licensed by Health Canada.